



# Labomat SPM

## Product information

### Product description

Labomat SPM is an enzymatic, low-foaming and economical instrument cleaner with high material compatibility, in concentrated form.

### Area of application

Labomat SPM is used for automated instrument reprocessing of surgical instruments and anaesthetic accessories.

### Instructions for use

Dosing takes place automatically via the dosing pump of the instrument dishwasher.

A temperature of 35-40°C is recommended for the inlet water. Temperatures significantly below 35°C can cause foaming. For the final rinsing, water (drinking water quality) is recommended.

The manufacturer's instructions for installation and maintenance must be observed. The parameters for the program sequence and the dosing settings must be adhered. Ion exchangers must be regenerated in a timely manner. For professional use only.

### Application concentration

0,2 - 1%

### Product range and specification

Item number	Container size	Container type	Packaging unit
600650	5 L	canister	3 x 5 L



# Labomat SPM

## Product information



### Medical device class

I

### Composition

- enzymes
- non-ionic surfactants
- glycols
- corrosion inhibitor

### Chemicophysical data

**Colour** slightly yellowish  
**Appearance** clear  
**State** liquid  
**pH-value (20°C)** 6,6 - 7,3  
**Density (20°C)** 1,025 - 1,036 g/ml

### Storage and shelf life

#### Shelf life

36 months

#### Shelf life after opening

12 months

Carefully reclose container after each use.

### Storage environment requirements

Store the product tightly closed in a well-ventilated location and protect it against unauthorised access. Observe the following storage temperature limits.

Description	Unit	Value
Temperature limits for short-term storage (min. - max.)	[°C]	0 - 40 <sup>1)</sup>
Ambient temperature for long-term storage according to shelf life (min. - max.)	[°C]	20 - 25

<sup>1)</sup> Short term storage (max. 24 h) up to max. 40 °C has no effect on the effectiveness of the product.

### Disposal

Disposal must be carried out in accordance with applicable national and local regulations. Please also note the information from the Material Safety Data Sheet.

### Information

Further information is available on the internet on our homepage:  
[www.drdeppe.de](http://www.drdeppe.de)

Special enquiries should be addressed to:  
[info@drdeppe.de](mailto:info@drdeppe.de)



**Laboratorium Dr. Deppe GmbH**  
 Hooghe Weg 35 · 47906 Kempen  
**T +49 (0) 2152 5565-0 W [www.drdeppe.de](http://www.drdeppe.de)**