



Product description

Spray Off wipes are ready-to-use, alcohol-based surface disinfectant wipes made of polypropylene material with a broad range of efficacy.

Area of application

For rapid disinfection of alcohol-resistant surfaces in the vicinity of patients, in hospitals, surgeries, outpatient clinics and care sectors as well as industry and public institutions. Excluded are surfaces of instruments or semi-critical and critical medical devices. As it is an alcohol-based disinfectant, it is not recommended for use on surfaces that are not resistant to alcohol, such as acrylic glass or polycarbonate.

Instructions for use

Observe the special hygiene requirements for the disinfection of surfaces that apply in the medical environment (in Germany: the recommendations of the KRINKO on the subject of surface disinfection, or the recommendations of the VAH for routine disinfection). The product is used as a wipe. It should be noted that the surface to be disinfected must be kept moist over the entire exposure time in order to take effect. The amount that has to be applied to the surface depends on the surface condition and ambient conditions such as temperature and humidity. In the food sector the areas with direct food contact must be rinsed thoroughly with water (drinking water quality) after disinfection. If necessary, wipe the surfaces dry with a hygienically clean wipe. For professional use only.

Exposure times

Spray Off wipes		Condition		Norm	Exposure time
		dirty	clean		
Bactericidal	with mechanical action	x		EN 13727*	1 min.
				EN 16615	
Yeasticidal		x		EN 13624*	
				EN 16615	
Fungicidal (C. albicans, A. brasiliensis)		x		EN 13624*	5 min.
Active against enveloped viruses (Vaccinia strain Elstree, BVDV strain NADL)			10% FBS	DVV/RKI*	1 min.
Limited spectrum of virucidal activity (Murine Norovirus, Adenovirus type 5)			10% FBS	DVV/RKI*	1 min.
Virucidal against Polyoma SV40			10% FBS	DVV/RKI*	5 min.
Virucidal (Murine Norovirus, Adenovirus type 5, Poliovirus type 1)			x	EN 14476*	5 min.

* tested on reference wipe

Product range and specification

Item number	Container size	Container type	Packaging unit
601351	150 wipes	can	12 x 150 wipes
601352	150 wipes	refill pack	12 x 150 wipes



601351



601352

Spray Off wipes

Product information



Medical device class

IIa

Biocide registration number

N-91346

Use biocides safely. Always read the label and product information before use.

Composition

The active ingredients contained in 100 g impregnating solution are:

63,50 g Ethanol

0,05 g Didecyltrimethylammonium chloride

Chemisch-physikalische Daten

Colour colourless

Appearance clear

State liquid

pH-value (20°C) 10,6 - 11,2

Density (20°C) 0,8829 - 0,8929 g/ml

Storage and shelf life

Shelf life

36 months

Shelf life after opening

28 days

Carefully reclose container after each use.

Only use damp wipes.

Storage environment requirements

Store the product tightly closed in a well-ventilated location and protect it against unauthorised access. Observe the following storage temperature limits.

Description	Unit	Value
Temperature limits for short-term storage (min. - max.)	[°C]	0 - 40 ¹⁾
Ambient temperature for long-term storage according to shelf life (min. - max.)	[°C]	20 - 25

¹⁾ Short term storage (max. 24 h) up to max. 40 °C has no effect on the effectiveness of the product.

Disposal

Disposal must be carried out in accordance with applicable national and local regulations. Please also note the information from the Material Safety Data Sheet.

Information

Further information is available on the internet on our homepage:

www.drdeppe.de

Special enquiries should be addressed to:

info@drdeppe.de



CE 0044



DrDeppe

Laboratorium Dr. Deppe GmbH
Hooghe Weg 35 · 47906 Kempen
T +49 (0) 2152 5565-0 **W** www.drdeppe.de