



Product description

Spray Off is a ready-to-use, liquid and alcohol-based surface disinfectant.

Area of application

It can be used for surface disinfection of medical devices and surfaces in a medical environment. Instrument surfaces or surfaces of semi-critical and critical medical devices are excluded. As it is an alcohol-based disinfectant, the surfaces must be alcohol resistant. Use on surfaces that are not alcohol resistant, such as acrylic glass, polycarbonate and aluminium, is not recommended without restrictions.

Instructions for use

The ready-to-use product can be sprayed, wiped on with a disposable cloth/mop, or in a combination of spraying/wiping. It should be noted that the surface to be disinfected must be kept moist over the entire exposure time in order to take effect. The amount of liquid that must be applied to the surface depends on the surface condition and ambient conditions such as temperature and humidity. For professional use only.

Exposure times

Spray Off		Condition		Norm	Exposure times
		dirty	clean		
Bactericidal	with mechanical action	x		EN 13727	1 min.
Yeasticidal		x		EN 16615*	
				EN 13624	
				EN 16615*	
Mycobactericidal (M. terrae, M. avium)		x		EN 14348	30 sec.
Fungicidal (C. albicans, A. brasiliensis)			x	EN 13624	5 min.
Active against enveloped viruses (Modified Vaccinia Ankara, BVDV)			10% FCS	DVV/RKI	1 min.
Limited spectrum of virucidal activity (Adenovirus Type 5, Murine Norovirus)			x	EN 14476	5 min.
Viruzidal (Poliovirus)			x	EN 14476	5 min.

* tested on reference wipe

Product range and specification

Item number	Container size	Container type	Packaging unit
600986	500 ml	bottle	18 x 500 ml
601330	1 L	bottle (spray head enclosed)	12 x 1 L
601334	5 L	canister	2 x 5 L
601332	10 L	canister	1 x 10 L



Spray Off

Product information

Medical device class

IIa

Biocide registration number

N-91345

Use biocides safely. Always read the label and product information before use.

Composition

The active ingredients contained in 100 g are:

63,50 g Ethanol

0,05 g Didecyldimethylammonium chloride

Chemicophysical data

Colour colourless

Appearance clear

State liquid

pH-value (20°C) 10,6 - 11,2

Density (20°C) 0,8829 - 0,8929 g/ml

Storage and shelf life

Shelf life

36 months

Shelf life after opening

12 months

Carefully reclose container after each use.

Storage environment requirements

Store the product tightly closed in a well-ventilated location and protect it against unauthorised access.

Observe the following storage temperature limits.

Description	Unit	Value
Temperature limits for short-term storage (min. - max.)	[°C]	0 - 40 ¹⁾
Ambient temperature for long-term storage according to shelf life (min. - max.)	[°C]	20 - 25

¹⁾ Short term storage (max. 24 h) up to max. 40 °C has no effect on the effectiveness of the product.

Disposal

Disposal must be carried out in accordance with applicable national and local regulations. Please also note the information from the Material Safety Data Sheet.

Information

Further information is available on the internet on our homepage:

www.drdeppe.de

Special enquiries should be addressed to:

info@drdeppe.de



CE 0044



Laboratorium Dr. Deppe GmbH
Hooghe Weg 35 · 47906 Kempen
T +49 (0) 2152 5565-0 **W** www.drdeppe.de