

Beta Guard Neu

Product information

Product description

Beta Guard Neu is a concentrated liquid, alcohol- and aldehyde-free surface disinfectant.

Area of application

For the disinfection of wipe-resistant surfaces in the patient environment, in clinics, practices, outpatient departments, and care areas, as well as in industry and public facilities. Excluded are surfaces of instruments, semi-critical and critical medical devices.

Spectrum of activity

- bactericidal (incl. MRSA), yeasticidal
- tuberculocidal, mycobactericidal
- active against enveloped viruses (incl. HBV, HIV, HCV, vaccinia-, influenza- and corona-viruses)
- limited spectrum of virucidal activity (incl. noro-, rota-, adeno-viruses)
- virucidal effective against Polyoma SV-40

Instructions for use

A diluted solution must be prepared before using the product. For this, the required amount of concentrate is mixed with water (drinking water quality). The prepared working solution should preferably be used by wiping with a disposable wipe. Spray disinfection should preferably be limited to areas that cannot be reached by wipe disinfection. When using the spray method, this should be used in a combination of spray/wipe method. The surface to be disinfected must be kept moist for the entire contact time in order to achieve effectiveness. For professional use only.



Product range and specification

item number	container size	container type	packaging unit	EAN/GTIN (each)	EAN/GTIN (packaging unit)
601487	2 L	canister	6 x 2 L	4260181252384	4260181252391
601050	5 L	canister	2 x 5 L	4260181254159	4260181251059
601052	10 L	canister	1 x 10 L	4260181254166	4260181251066

Beta Guard Neu

Product information

Medical device class

IIa

Biocide registration number

N-77462

Use disinfectant safely. Always read the label and product information before use.

Composition

The active ingredients contained in 100 g:

- 7.0 g Alkyl (C12-16) dimethylbenzyl ammonium chloride
- 5.0 g N-(3-Aminopropyl)-N-dodecylpropane-1,3-diamine
- 2.0 g 2-Phenoxyethanol

Storage and shelf life

Shelf life

36 months

Shelf life after opening

12 months

Close the container carefully after each use.


Store the product in a well-ventilated place in a sealed condition and protect it from unauthorized access.

Exposure times

Beta Guard Neu		condition		norm	concentration	exposure time (min.)
		dirty	clean			
bactericidal	with mechanical action	x		EN 13727	0,5 %	5
				EN 16615*		
yeastical		x		EN 13624	1,5 %	5
				EN 16615*		
bactericidal	without mechanical action	x		EN 13727 / EN 1276	1,5 %	5
				EN 13697		
yeastical		x		EN 13624 / EN 1650	1,5 %	5
				EN 13697		
tuberculocidal/ mycobactericidal (M. terrae, M. avium)		x		EN 14348	2 % 1 %	30** 60**
active against enveloped viruses (Modified Vaccinia Ankara)		x		EN 14476	3 %	5
limited spectrum of virucidal activity:		x		EN 14476	4 %	60**
					4 %	60**
Adenovirus type 5		x		EN 14476	3 %	60**
Murine Norovirus		x			1 %	30**
Rotavirus strain Wa			without		1 %	30**
virucidal effective against: Polyomavirus SV-40		x		EN 14476	2 %	60**

* tested on reference wipe ** not for near patient surfaces

Information about this product



For further information, please visit our website www.drdeppe.de

For specific inquiries, please contact info@drdeppe.de

CE 0044



Laboratorium Dr. Deppe GmbH
Hooghe Weg 35 · 47906 Kempen
T +49 (0) 2152 5565-0 W www.drdeppe.de