

Beta Guard rfu wipes

Product information

Product description

Beta Guard rfu wipes are ready-to-use surface disinfectant wipes that are alcohol and aldehyde-free. The used wipe material is polypropylene.

Area of application

For the disinfection of wipe-resistant device surfaces and objects in the patient environment, as well as surfaces in industrial and food areas. Excluded are surfaces of instruments, semi-critical and critical medical devices.

Spectrum of activity

- bactericidal (incl. MRSA), yeasticidal
- tuberculocidal
- active against enveloped viruses (incl. HBV, HIV, HCV, vaccinia-, influenza- and corona-viruses)
- limited spectrum of virucidal activity (incl. noro-, rota-, adeno-viruses)

Instructions for use

The ready-to-use product is used by wiping. It is essential to keep the surface to be disinfected moist throughout the entire exposure times to achieve effectiveness.



Product range and specification

item number	container size	wipe size	container type	packaging unit	EAN/GTIN (each)	EAN/GTIN (packaging unit)
601111	150 wipes	15 x 18 cm	can	12 x 150 wipes	4260181253992	4260181251110
601113	150 wipes	15 x 18 cm	can-refill pack	12 x 150 wipes	4260181254005	4260181251127

Beta Guard rfu wipes

Product information

Biocide registration number

N-83979

Use disinfectant safely. Always read the label and product information before use.

Composition

The active ingredients contained in 100 g soaking solution:

0.1 g 2-Phenoxyethanol

0.35 g Alkyl(C12-16) dimethyl benzyl ammonium chloride

0.25 g N-(3-Aminopropyl)-N-dodecylpropane-1,3-diamine

Storage and shelf life

Shelf life

can and can refillpack 36 months

Shelf life after opening

28 days

Close the container carefully after each use.

Only use moist wipes.

Store the product in a well-ventilated place in a sealed condition and protect it from unauthorized access.

Exposure times

Beta Guard rfu wipes	condition		norm	Exposure time
	dirty	clean		
bactericidal	x		EN 13727* EN 16615	1 min.
yeastcidal	x		EN 13624* EN 16615	
tuberculocidal (M. terrae)	x		EN 14348*	60 min.**
active against enveloped viruses	x		EN 14476*	1 min.
Modified Vaccinia Ankara		x	mod. EN 16615	
limited spectrum of virucidal activity (Adenovirus type 5, Murine Norovirus)	x		EN 14476*	60 min.**

* tested on soaking solution

** not for near patient surfaces

Informations

about this product



For further information, please visit our website
www.drdeppe.de

For specific inquiries, please contact
info@drdeppe.de



Laboratorium Dr. Deppe GmbH
Hooghe Weg 35 · 47906 Kempen
T +49 (0) 2152 5565-0 W www.drdeppe.de