


Product information

Instru Extra is an aldehyde- and phenol-free instrument disinfectant concentrate for the reprocessing of surgical instruments, thermolabile anaesthesia accessories, and rigid endoscopes in immersion or ultrasonic bath. Also suitable for drills or rotating instruments.

For all areas where medical instruments are reprocessed using immersion or ultrasonic bath methods, Instru Extra is not suitable for the final disinfection of semi-critical instruments.



- bactericidal (including MRSA), yeasticidal
- tuberculocidal
- active against enveloped viruses (including HBV, HIV, HCV, vaccinia-, influenza-, and corona viruses)
- instrument disinfection in an ultrasonic bath

Insert instruments after cleaning. All surfaces and cavities must be covered or filled. Please follow the instructions for application concentration and contact times. Remove instruments, rinse thoroughly with water, and dry. Follow the hygiene requirements for the reprocessing of medical devices (in Germany: RKI, KRINKO, MPBetriebV). Follow the instructions of the instrument and equipment manufacturers.

Similar to the immersion method. Please follow the contact times specified for the ultrasonic application. Follow the instructions of the instrument and equipment manufacturers.

Instru Extra is not compatible with aldehyde-based disinfectants. For initial use and product change, thoroughly clean and rinse the instruments and accessories, trays, reprocessing equipment, and application aids. Instru Extra working solution is resistant to protein contamination and can be used multiple times if necessary. If cloudiness or visible contamination occurs, it is recommended to replace the working solution. For professional use only.

p. 2 ►

Instru Extra

Product information

Medical device class

IIb

Composition

The active ingredients contained in 100 g:

5 g Alkyldimethylethylammoniummethylsulfate

1 g Polyhexamethylenbiguanide hydrochloride

1.8 g Cocosporylenguanidiniumdiacetate

Storage and shelf life

Shelf life

36 months

Shelf life after opening

12 months

Close the container carefully after each use.

Store the product in a well-ventilated place in a sealed condition and protect it from unauthorized access.

Exposure times

Instru Extra	condition		Norm	concentration	exposure time
	dirty	clean			
bactericidal		x	EN 13727 EN 14561 / VAH 15	3%	15 min.
yeastcidal		x	EN 13624 EN 14562 / VAH 15		
tuberculocidal (M. terrae)		x	EN 14348 EN 14563		
active against enveloped viruses acc. to DVV/RKI (Vaccinia strain Elstree, BVDV)	10% FBS		DVV/RKI	1%	15 min.
active against enveloped viruses acc. to EN 14476 (Modified Vaccinia Ankara)		x	EN 14476		
instrument disinfection in an ultrasonic bath	x		VAH 15	3%	5 min.

Information

about this product



For further information,
please visit our website

www.drdeppe.de

For specific inquiries,
please contact

info@drdeppe.de

CE 0044



Laboratorium Dr. Deppe GmbH

Hooghe Weg 35 · 47906 Kempen

T +49 (0) 2152 5565-0 W www.drdeppe.de