

Labomat KS

Product information

Product description

Labomat KS is a ready-to-use, phosphate- and surfactant-free rinsing agent for disinfection dishwashers.

Area of application

Labomat KS is used to neutralise alkaline substances on medical instruments and laboratory equipment after fully automated cleaning.

Application concentration

Setting of the dosing pump: 0.1-0.5%

Instructions for use

Dosing takes place automatically via the dosing pump in the dishwasher. Please follow the instructions of the relevant device manufacturer. Do not mix with other products. When changing products, rinse the dosing system including the suction tubing with water. Labomat KS should not be used on objects made of anodised aluminium, chrome-plated or nickel-plated parts. For professional use only.



Product range and specification

item number	container size	container type	packaging unit	EAN/GTIN (each)	EAN/GTIN (packaging unit)
600021	5 L	canister	3 x 5 L	4260181253961	4260181251387

**DrDeppe**

Labomat KS

Product information

Product type

Detergent

Composition

Ingredients

Citric acid

Storage and shelf life

Shelf life

36 months

Shelf life after opening

12 months

Close the container carefully after each use.

Store the product in a well-ventilated place in a sealed condition and protect it from unauthorized access.

Information

about this product

For further information,
please visit our website**www.drdeppe.de**For specific inquiries,
please contact**info@drdeppe.de****DrDeppe****Laboratorium Dr. Deppe GmbH**

Hooghe Weg 35 · 47906 Kempen

T +49 (0) 2152 5565-0 **W** www.drdeppe.de