

# Spray In GA

## Product information

## Product description

Spray In GA is a ready-to-use, liquid, low-alcoholic disinfectant for surface disinfection. The product shows very good material compatibility under the specified exposure times. Spray In GA is available in the variants Neutral (without fragrance), Melone (fruity fragrance) and Aktiv (fresh fragrance).

## Area of application

For rapid disinfection of alcohol resistant surfaces in the vicinity of patients, in clinics, surgeries, outpatient departments and care areas, as well as in industry and public facilities.



## Spectrum of activity

- bactericidal (incl. MRSA), yeasticidal
- tuberculocidal, mycobactericidal
- active against enveloped viruses (incl. HBV, HIV, HCV, vaccinia-, influenza- and corona-viruses)
- limited spectrum of virucidal activity (incl. noro-, rota-, adeno-viruses)
- effective against Polyoma SV-40

## Instructions for use

The ready-to-use product is preferably to be used in the wiping method using a disposable wipe. Spray disinfection should preferably be limited to areas that are not reachable by wiping disinfection. When using the spray method, it should be combined with wiping method. It is important to keep the surface to be disinfected moist throughout the entire contact time to achieve effectiveness.

## Product range and specification

item number	container size	container type	packaging unit	EAN/GTIN (each)	EAN/GTIN (packaging unit)
601611 (Neutral)	500 ml	bottle	18 x 500 ml	4260181253459	4260181253466
601612 (Melone)	500 ml	bottle	18 x 500 ml	4260181253473	4260181253480
601613 (Aktiv)	500 ml	bottle	18 x 500 ml	4260181253497	4260181253503

# Spray In GA

## Product information

### Biocide registration number

N-109576

Use disinfectant safely. Always read the label and product information before use.

### Composition

**The active ingredients contained in 100 g:**

20 g Propan-1-ol

0.2 g Alkylbenzyltrimethylammonium chloride

### Storage and shelf life

#### Shelf life

36 months

#### Shelf life after opening

12 months

Carefully reclose container after each use.

Store the product tightly closed in a well-ventilated location and protect it against unauthorised access.

### Exposure times

Spray In GA		condition		norm	exposure time
		dirty	clean		
bactericidal	with mechanical action	x		EN 13727 EN 16615*	1 min.
yeastcidal		x		EN 13624 EN 16615*	
mycobactericidal (M. terrae, M. avium)		x		VAH 9.1 VAH 14.1*	2 min.
active against enveloped viruses (Modified Vaccinia Ankara)			x	EN 14476 modified	1 min.
limited spectrum of virucidal activity:			x	EN 14476	30 min.**
Adenovirus type 5			x	EN 14476	30 min.**
Murine Norovirus			x	EN 14476	1 min.
Rotavirus strain WA			with-out	EN 14476	30 sec.
virucidal effective against Polyoma SV-40			x	EN 14476	5 min.

\* tested on reference wipe

\*\* not for near patient surfaces

### Information

about this product



For further information, please visit our website

**[www.drdeppe.de](http://www.drdeppe.de)**

For specific inquiries, please contact

**[info@drdeppe.de](mailto:info@drdeppe.de)**



**Laboratorium Dr. Deppe GmbH**

Hooghe Weg 35 · 47906 Kempen

**T** +49 (0) 2152 5565-0 **W** [www.drdeppe.de](http://www.drdeppe.de)