

# Spray In Neu

## Product information

## Product description

Spray In Neu is a ready-to-use liquid, alcoholic and aldehyde-free surface disinfectant.

## Area of application

For rapid disinfection of alcohol-resistant surfaces in patient environments, clinics, practices, outpatient departments, and care areas, as well as in industry and public facilities. Excluded are surfaces of instruments, semi-critical and critical medical devices.



## Spectrum of activity

- bactericidal (incl. MRSA), yeasticidal
- tuberculocidal
- active against enveloped viruses (incl. HBV, HIV, HCV, vaccinia-, influenza- and corona-viruses)
- limited spectrum of virucidal activity (incl. noro-, rota-, adeno-viruses)

## Instructions for use

The ready-to-use product is preferably to be used in the wiping method using a disposable wipe. Spray disinfection should preferably be limited to areas that are not reachable by wiping disinfection. When using the spray method, it should be combined with wiping method. It is important to keep the surface to be disinfected moist throughout the entire contact time to achieve effectiveness.

## Product range and specification

item number	container size	container type	packaging unit	EAN/GTIN (each)	EAN/GTIN (packaging unit)
601450	1 L	club bottle	12 x 1 L	4260181254234	4260181250991
601454	5 L	canister	2 x 5 L	4260181254241	4260181251004
601452	10 L	canister	1 x 10 L	4260181254258	4260181251011

# Spray In Neu

Product information

## Medical device class

IIa

## Biocide registration number

N-60215

Use disinfectant safely. Always read the label and product information before use.

## Composition

**The active ingredients contained in 100 g:**

0.05 g Didecyltrimethylammonium chloride

43 g Ethanol

## Storage and shelf life

**Shelf life**

36 months

**Shelf life after opening**

12 months

Close the container carefully after each use.

Store the product in a well-ventilated place in a sealed condition and protect it from unauthorized access.

## Exposure times

Spray In Neu		condition		norm	exposure time
		dirty	clean		
bactericidal	with mechanical action	x		EN 13727 EN 16615*	5 min.
yeastcidal		x		EN 13624 EN 16615*	
bactericidal	without mechanical action	x		EN 13727 EN 13697	
yeastcidal		x		EN 13624 EN 13697	
tuberculocidal (M. terrae)		x		EN 14348	1 min.
active against enveloped viruses (Vaccinia strain Elstree, BVDV strain NADL)		10% FBS		DVV/RKI	30 sec.
active against enveloped viruses (Modified Vaccinia Ankara)			x	EN 14476	1 min.
limited spectrum of virucidal activity:			x	EN 14476	15 min.**
Adenovirus type 5			x	EN 14476	15 min.**
Murine Norovirus			x	EN 14476	5 min.
Rotavirus strain Wa			x	EN 14476	5 min.

\* tested on reference wipe

\*\* not for near patient surfaces

## Information

about this product



For further information, please visit our website

**[www.drdeppe.de](http://www.drdeppe.de)**

For specific inquiries, please contact

**[info@drdeppe.de](mailto:info@drdeppe.de)**

CE 0044



**Laboratorium Dr. Deppe GmbH**

Hooghe Weg 35 · 47906 Kempen

T +49 (0) 2152 5565-0 W [www.drdeppe.de](http://www.drdeppe.de)