

Beta Guard rfu

Product information

Product description

Beta Guard rfu is a ready-to-use, liquid, alcohol- and aldehyde-free surface disinfectant.

Area of application

It is suitable for rapid disinfection of wipe-resistant surfaces in patient environments, clinics, practices, outpatient departments, and care areas, as well as in industry and public facilities.

Spectrum of activity

- bactericidal (incl. MRSA), yeasticidal
- tuberculocidal
- active against enveloped viruses (incl. HBV, HIV, HCV, vaccinia-, influenza- and corona-viruses)
- limited spectrum of virucidal activity (incl. noro-, rota-, adeno-viruses)

Instructions for use

The ready-to-use product is preferably to be used in the wiping method using a disposable wipe. Spray disinfection should preferably be limited to areas that are not reachable by wiping disinfection. When using the spray method, it should be combined with wiping method. It is important to keep the surface to be disinfected moist throughout the entire contact time to achieve effectiveness.



Product range and specification

item number	container size	container type	packaging unit	PCN (each)	EAN/GTIN (each)	EAN/GTIN (packaging unit)
601544*	500 ml	bottle	18 × 500 ml	18312289	4260181253060	4260181253077
601110	1 L	club bottle	12 × 1 L	18312929	4260181254180	4260181251028
601112	5 L	canister	2 × 5 L	18312941	4260181254197	4260181251035
601114	10 L	canister	1 × 10 L	18312958	4260181254203	4260181251042

* item is a biocidal product

Beta Guard rfu

Product information

Medical device class

IIa (except for item no. 601544 - biocidal product only)

Biocide registration number

N-83978

Use disinfectant safely. Always read the label and product information before use.

Composition

The active ingredients contained in 100 g:

0.1 g 2-Phenoxyethanol

0.35 g Alkyl (C12-16) dimethylbenzyl ammonium chloride

0.25 g N-(3-Aminopropyl)-N-dodecylpropane-1,3-diamine

Storage and shelf life

Shelf life:

36 months

Shelf life after opening:

12 months

Close the container carefully after each use.

Store the product in a well-ventilated place in a sealed condition and protect it from unauthorized access.

Exposure times

Beta Guard rfu spectrum of activity	condition (dirty/clean)	exposure time	norm
bactericidal	with mechanical action dirty	1 minute	EN 13727 EN 16615*
yeasticidal	with mechanical action dirty	1 minute	EN 13624 EN 16615*
bactericidal	without mechanical action dirty	1 minute	EN 13727 EN 1276*** EN 13697
yeasticidal	without mechanical action dirty	1 minute	EN 13624 EN 1650*** EN 13697
tuberculocidal (M. terrae)	dirty	60 minutes**	EN 14348
active against enveloped viruses (Modified Vaccinia Ankara)	dirty	1 minute	EN 14476, virus-4-field-test****
limited spectrum of virucidal activity (Adenovirus type 5, Murine Norovirus)	dirty	60 minutes**	EN 14476


* tested on reference wipe

** not for near patient surfaces

*** food sector

**** according to EN 16615

Informationen about this product



For further information, please visit our website **www.drdeppe.de**

For specific inquiries, please contact **info@drdeppe.de**