

# Endo Star Neu

## Product information

### Product description

Endo Star Neu is an aldehyde- and phenol-free instrument disinfectant concentrate for the reprocessing of surgical instruments, as well as thermolabile anaesthesia equipment and rigid endoscopes in immersion or ultrasonic bath. Also suitable for drills or rotating instruments.

### Area of application

For all areas where medical instruments are reprocessed using immersion or ultrasonic bath methods. Endo Star Neu is not suitable for the final disinfection of semicritical instruments. For professional use only.



### Spectrum of activity

- bactericidal (including MRSA), yeasticidal
- tuberculocidal
- active against enveloped viruses (including HBV, HIV, HCV, vaccinia-, influenza-, and corona viruses)
- instrument disinfection in an ultrasonic bath

### Instructions for use

#### Immersion method

Insert instruments after cleaning. All surfaces and cavities must be covered or filled. Please follow the application concentration and contact times. Remove instruments, rinse thoroughly with water (potable quality), and dry. Follow the hygiene requirements for the reprocessing of medical devices (in Germany: RKI, KRINKO, MPBetreibV). Follow the instructions of the instrument and equipment manufacturers.

#### Ultrasound application

Similar to the immersion method. Adhere to the contact times specified for ultrasound application.

### Product change

Endo Star Neu is not compatible with aldehyde-based disinfectants. For initial use and product change, thoroughly clean and rinse the instruments and accessories, trays, reprocessing equipment, and application aids. The working solution is resistant to protein and can potentially be reused multiple times. If cloudiness or visible contamination occurs, it is recommended to replace the working solution.

### Product range and specification

item number	container size	container type	packaging unit	PCN (each)	EAN/GTIN (each)	EAN/GTIN (packaging unit)
601486	2 L	canister	6 × 2 L	18313231	4260181252360	4260181252377
600080	5 L	canister	2 × 5 L	18313225	4260181253800	4260181251233

# Endo Star Neu

Product information

## Medical device class

IIb

## Composition

**The active ingredients contained in 100 g:**

5 g Alkyldimethylethylammoniummethylsulfate

1 g Polyhexamethylenebiguanide hydrochloride

1.8 g Cocosporylenguanidiniumdiacetate

## Storage and shelf life

**Shelf life:**

36 months

**Shelf life after opening:**

12 months

Close the container carefully after each use.

Store the product in a well-ventilated place in a sealed condition and protect it from unauthorized access.

## Exposure times

Endo Star Neu spectrum of activity	condition (dirty/clean)	concentration	exposure time	norm
bactericidal	clean	3%	15 min.	EN 13727 EN 14561 (VAH 15)
yeastcidal	clean			EN 13624 EN 14562 (VAH 15)
tuberculocidal (M. terrae)	clean			EN 14348 EN 14563
active against enveloped viruses (Vaccinia strain Elstree, BVDV)	10% FBS			DVV/RKI
active against enveloped viruses (Modified Vaccinia Ankara)	clean	1%	15 min.	EN 14476
instrument disinfection in ultrasonic bath	dirty	3%	5 min.	VAH 15

## Information about this product



For further information, please visit our website [www.drdeppe.de](http://www.drdeppe.de)

For specific inquiries, please contact [info@drdeppe.de](mailto:info@drdeppe.de)

CE 0044



**Laboratorium Dr. Deppe GmbH**  
Hooghe Weg 35 · 47906 Kempen  
T +49 (0) 2152 5565-0 W [www.drdeppe.de](http://www.drdeppe.de)