

# Labomat SPM

## Product information

### Product description

Labomat SPM is a concentrated enzymatic, low-foaming and economical instrument cleaner with high material compatibility.

### Area of application

Labomat SPM is used for the mechanical reprocessing of surgical instruments and anaesthesia accessories. The enzymes used in Labomat SPM are temperaturestable up to 70°C.

### Application concentration

0.2-1%

### Instructions for use

The dosage is carried out using the dosing pumps in the instrument dishwasher. A temperature of 35-40°C is recommended for the inlet water. Temperatures significantly below 35°C may cause foaming. The instructions of the device manufacturers regarding installation, operation, and maintenance must be observed, especially regarding the quality requirements of the process water being used. The validated parameters for the program sequence and dosing settings must be adhered to. Ion exchangers must be regenerated in a timely manner. For professional use only.



### Product range and specification

item number	container size	container type	packaging unit	PCN (each)	EAN/GTIN (each)	EAN/GTIN (packaging unit)
600650	5 L	flat-top canister	3 × 5 L	18313478	4260181253671	4260181251394

# Labomat SPM

Product information

## Medical Device Class

I

## Composition

### Ingredients

Enzymes, non-ionic surfactants, glycols, purified water

## Storage and shelf life

### Shelf life:

36 months

### Shelf life after opening:

12 months

Close the container carefully after each use.

Store the product in a well-ventilated place in a sealed condition and protect it from unauthorized access.

### Information about this product



For further  
information, please  
visit our website  
**[www.drdeppe.de](http://www.drdeppe.de)**

For specific inquiries,  
please contact  
**[info@drdeppe.de](mailto:info@drdeppe.de)**