

# Spray Off wipes

Product information

## Product description

Spray Off wipes are ready-to-use, alcohol-based, aldehyde-free surface disinfectant wipes. The used wipe material is polypropylene.

## Area of application

For surface disinfection of alcohol-resistant surfaces in the patient-near environment, clinics, practices, outpatient clinics and care facilities, as well as in industry and public facilities. Surfaces of instruments resp. semi-critical and critical medical devices are excluded.

## Spectrum of activity

- bactericidal (incl. MRSA), yeasticidal
- tuberculocidal, mycobactericidal
- fungicidal
- active against enveloped viruses (incl. HBV, HIV, HCV, vaccinia-, influenza- and corona-viruses)
- limited spectrum of virucidal activity (incl. noro-, rota-, adeno-viruses)
- virucidal effective against Polyoma SV-40
- virucidal

## Instructions for use

The ready-to-use product is used by wiping. It is essential to keep the surface to be disinfected moist throughout the entire exposure times to achieve effectiveness.



## Product range and specification

item number	container size	wipe size	container type	packaging unit	PCN (each)	EAN/GTIN (each)	EAN/GTIN (packaging unit)
601351	150 wipes	15 × 18 cm	can	12 × 150 wipes	18313202	4260181254425	4260181254524
601352	150 wipes	15 × 18 cm	can-refill pack	12 × 150 wipes	18313219	4260181254432	4260181254531



# Spray Off wipes

Product information

## Medical device class

IIa

## Biocide registration number

N-91346

Use disinfectant safely. Always read the label and product information before use.

## Composition

The active ingredients contained in 100 g soaking solution:

63.50 g Ethanol

0.05 g Didecyldimethylammonium chloride

## Storage and shelf life

### Shelf life:

36 months

### Shelf life after opening:

28 days

Close the container carefully after each use. Only use moist wipes.

Store the product in a well-ventilated place in a sealed condition and protect it from unauthorized access.

## Exposure times

Spray Off wipes (Polypropylene) spectrum of activity	condition (dirty/clean)	exposure time	norm
bactericidal	with mechanical action dirty	1 minute	EN 13727* EN 16615
yeasticidal	with mechanical action dirty	1 minute	EN 13624* EN 16615
tuberculocidal (M. terrae), mycobactericidal (M. terrae, M. avium)	dirty	1 minute	EN 14348*
fungicidal (C. albicans, A. brasiliensis)	dirty	5 minutes	EN 13624*
active against enveloped viruses (Vaccinia strain Elstree, BVDV strain NADL)	10% FBS	1 minute	DVV/RKI*
active against enveloped viruses (Modified Vaccinia Ankara)	clean	1 minute	EN 14476*
limited spectrum of virucidal activity (Murine Norovirus, Adenovirus type 5)	10% FBS	1 minute	DVV/RKI*
virucidal effective against Polyoma SV-40	10% FBS	5 minutes	DVV/RKI*
virucidal (Murine Norovirus, Adenovirus type 5, Poliovirus type 1)	clean	5 minutes	EN 14476*

\*tested on soaking solution

## Information about this product



For further information, please visit our website [www.drdeppe.de](http://www.drdeppe.de)

For specific inquiries, please contact [info@drdeppe.de](mailto:info@drdeppe.de)

CE 0044



Laboratorium Dr. Deppe GmbH  
Hooghe Weg 35 · 47906 Kempen

T +49 (0) 2152 5565-0 W [www.drdeppe.de](http://www.drdeppe.de)