



Spray In QF wipes (cellulose)

Surface disinfection

EN-Instructions for use

Thank you for purchasing this product from Laboratorium Dr. Deppe GmbH.

If you have any questions about the contents of these instructions for use or the application of the product, please contact us. Your Laboratorium Dr. Deppe Team.



Laboratorium Dr. Deppe GmbH
Hooghe Weg 35
D - 47906 Kempen
+49 (0) 2152 5565-0
kontakt@drdeppe.de
www.drdeppe.de

CE 0044

Instructions for use version: 2.0_EN
Date of issue: 05.03.2026

REF 480079-1

About these instructions for use

These instructions for use contain the necessary information for the correct use of Spray In QF wipes (referred to below as the "product"). The instructions for use must be kept in the immediate vicinity of the product and must be available to the user at all times. Anyone using the product must have read and understood the instructions for use in full, in particular the section on safety and the accompanying documents. Report all serious incidents (damage, injuries, infections, etc.) that occur in connection with the product to the manufacturer (see page 2) and the competent authority of the Member State in which the user and/or patient is established.

Accompanying documents

The following documents must be complied with in addition to the instructions for use: Labels on the product packaging, product data sheets, material safety data sheets, operating instructions, information on the shelf life of open containers. The documents can be downloaded from the Laboratorium Dr. Deppe GmbH website or requested by e-mail.

Explanation of symbols

These instructions for use and the labelling on the product use symbols to draw your attention to certain information and hazards.

Symbol	Description
	Manufacturer
	Use-by date
	Lot number
	Catalogue number, DrDeppe item number
	Consult instructions for use
	Medical device
	Unique Device Identifier
	The letters CE stand for European Conformity; this symbolises the compliance of the product with the applicable requirements that the European Community places on the manufacturer.
	Flammable

Safety

The following section provides information on the safe use of the product. Read and follow all safety instructions and contact the manufacturer if anything is unclear.

Personal protective equipment*

Eye protection: Safety goggles are recommended when decanting. Not required for normal handling.

Skin and body protection: Wear suitable protective clothing when working.

Hand protection: Not required for normal handling.

Respiratory protection: Wear suitable respiratory protective equipment if ventilation is inadequate. Not required for normal handling. Short-term exposure. A-P2.

*Section 8.2. Exposure controls / Safety data sheet

Safety information

The following safety instructions apply when handling the product

Do not use the product on humans. Do not use the product for the disinfection of instruments, semi-critical and critical medical devices. Check that the packaging is intact before use. Do not use the product if the packaging is damaged. Do not use the product after the expiry date. Wear the necessary personal protective equipment. Use only fully saturated disinfectant wipes. Be aware of the risk of explosion and fire when using alcohol-based disinfectants on large surfaces. Ensure good ventilation during and after use. Do not eat, drink or smoke when using the device. Protect the product from direct sunlight and keep it away from sources of ignition, food, drink and animal feed. Avoid contact with the skin, eyes and mucous membranes. Do not inhale mists, vapours or aerosols. In case of accidental release, carefully clean and ventilate the contaminated area. Avoid releasing it into the environment. Wash your hands carefully after use. Remove clothing that has come into contact with the product immediately after use and wash it before wearing it again. Repeated contact with the product can lead to chapped or cracked skin. Take precautionary measures against static discharge. Keep out of reach of children.

Safety measures

Take the following measures in the event of improper use or accidents with the product.

First-aid measures

General: If you are unsure of anything or if symptoms persist, consult a doctor. After inhalation: Move the person to fresh air and ensure unobstructed breathing.

After skin contact: The product is not considered to be irritating to the skin.

Acute and delayed symptoms: Repeated exposure may cause skin dryness or cracking.

After eye contact: Rinse eyes with water as a precaution.

After ingestion: Call a poison centre or doctor if you feel unwell.

Firefighting measures

Comply with the following instructions and measures in the event of fire. Suitable extinguishing media: Suitable extinguishing media for surrounding fires (water spray jet, dry extinguishing powder, foam, carbon dioxide). Unsuitable extinguishing media: Full jet of water.

Special hazards arising from the substance or mixture

In the event of fire, toxic fumes and carbon oxides (CO, CO₂) may be released as hazardous decomposition products.

Advice for firefighters

Protection during firefighting: Do not attempt to work without suitable protective equipment. Use a self-contained breathing apparatus and full protective clothing.

Disposal of extinguishing media: Do not allow the extinguishing water to enter the sewage system or watercourses. The disposal of extinguishing water must be carried out in accordance with official regulations.

Accidental release measures

Comply with the following instructions and measures in the event of accidental release of the product.

Personal safety precautions

For non-emergency personnel: Ventilate the contaminated area. void open flames, sparks and smoking.

For emergency responders: Do not attempt to work without suitable protective equipment.

Environmental precautions

If the product enters the sewage system or public waters, the authorities must be notified. Avoid penetration into the subsoil. Prevent the product from entering the sewage system and public waters.

Methods and material for containment and cleaning up

For containment: Collect spillage.

Cleaning procedure: Absorb the spilled liquid with absorbent material. Pick up the absorbent material using mechanical methods (wiping up, sweeping up) and collect it in suitable containers for disposal.

Disposal: Disposal must be carried out in accordance with official regulations.

Intended purpose

The product is an alcohol-based surface disinfectant for disinfection in particular of surfaces of non-invasive medical devices and surfaces in medical environments.

Requirements for the user

The product is intended exclusively for use by trained medical personnel and cleaning specialists from the healthcare sector. Anyone using the product must have read and understood the entire instructions for use and in particular the section on safety before use.

Requirements for the usage environment

The product is intended for use within the following temperature ranges.

Outside the temperature limits there are no data on efficacy.

Recommended ambient temperature for use (min. - max.): 20 - 25°C

Material compatibility

The product shows good material compatibility under the specified exposure times. For surfaces that are not resistant to alcohol, such as acrylic glass and polycarbonate, optical changes such as faint traces of haze or a matt sheen may occur. We recommend testing material compatibility for any specific application on an inconspicuous area before the first application.

Product properties

The product can be used as a disinfectant for professional surface disinfection of non-invasive medical devices and surfaces in medical environments. It serves to prevent diseases caused by viruses and bacteria by reducing the number of pathogenic germs and thus minimising the likelihood of infection.

Composition

The following information refers to 100 g of the product:

Ethanol: 50,00 g

Chemical and physical properties

Colour: clear

State: liquid

pH-value: 7 - 9

Density: 0,91 - 0,9250 g/ml

Contact Times

Product effectiveness is guaranteed in compliance with the specified contact times and temperature limits. The amount that needs to be applied to the surface depends on the application conditions.

Influencing factors: Increase the exposure time for rough surfaces and low temperatures. Keep the surface moistened with working solution for the entire exposure time for a safe effect.

Efficacy: bactericidal (EN 13727*; EN 1276*EN 16615)/yeasticidal (EN 13624*; EN 1650*; EN 16615) clean condition 1 min.; active against enveloped viruses (Modified Vaccinia Ankara/EN 14476*; virus-4-field-test**) clean condition 1 min.; limited spectrum of virucidal activity (Adenovirus type 5, Murine Norovirus/ EN 14476) clean condition 1 min.; tuberculocidal (m. terrae/EN 14348*) clean condition 1 min.; mycobactericidal (M. terrae, M. avium/EN 14348*), clean condition 5 min.

* tested on reference wipe ** based on EN 16615

Use

Dosing the product correctly

Carry out the following steps to disinfect surfaces with the product. Remove visible soiling from the surface to be disinfected with a soft, clean wipe.

Wipe disinfection with ready-to-use disinfectant wipe

Open the closure flap, remove a disinfectant wipe from the dispenser, and securely close the flap again. The disinfectant wipes must be fully saturated with the soaking solution. Test material compatibility on an inconspicuous area, taking into account the planned contact time. Carefully wipe the entire surface with the disinfectant wipe until it is completely wetted with the ready-to-use solution/soaking solution. If necessary, wipe the surface again with a fresh disinfectant wipe. Dispose of the disinfectant wipe properly.

Storage and disposal

Requirements for the storage environment

Store the product in a well-ventilated place in a sealed condition and protect it from unauthorised access.

Shelf life

Do not use the product after the expiry date indicated on the packaging and dispose of it properly. Ready-to-use disinfectant wipes can be used for up to 28 days after opening if kept sealed.

Disposal

Disposal must be carried out in accordance with the applicable national and local regulations. The instructions in the safety data sheet for the disinfectant wipes must also be complied with.