

Spray Off

Product information

Product description

Spray Off is a liquid ready-to-use, alcoholic and aldehyde-free surface disinfectant.

Area of application

For the disinfection of alcohol-resistant surfaces in the patient environment, in clinics, practices, outpatient departments, and care areas, as well as in industry and public facilities. Excluded are surfaces of instruments, semi-critical and critical medical devices.



Spectrum of activity

- bactericidal (incl. MRSA), yeasticidal
- fungicidal
- active against enveloped viruses (incl. HBV, HIV, HCV, vaccinia-, influenza- and corona-viruses)
- limited spectrum of virucidal activity (incl. noro-, rota-, adeno-viruses)
- virucidal effective against Polyoma SV-40
- virucidal

Instructions for use

The ready-to-use product is preferably to be used in the wiping method using a disposable wipe. Spray disinfection should preferably be limited to areas that are not reachable by wiping disinfection. When using the spray method, it should be combined with wiping method. It is important to keep the surface to be disinfected moist throughout the entire contact time to achieve effectiveness.

Product range and specification

item number	container size	container type	packaging unit	PCN (each)	EAN/GTIN (each)	EAN/GTIN (packaging unit)
600986	500 ml	bottle	18 x 500 ml	18313136	4260181253794	4260181250762
601330	1 L	club bottle	12 x 1 L	18313142	4260181254357	4260181250779
601334	5 L	canister	2 x 5 L	18313159	4260181254371	4260181250786
601332	10 L	canister	1 x 10 L	18313165	4260181254388	4260181250793

Spray Off

Product information

Medical device class

IIa

Biocide registration number

N-91345

Use disinfectant safely. Always read the label and product information before use.

Composition

The active ingredients contained in 100 g:

63.50 g Ethanol

0.05 g Didecyldimethylammonium chloride

Storage and shelf life

Shelf life:

36 months

Shelf life after opening:

12 months

Carefully reclose container after each use.

Store the product tightly closed in a well-ventilated location and protect it against unauthorised access.

Exposure times

Spray Off spectrum of activity	condition (dirty/clean)	exposure time	norm
bactericidal	with mechanical action dirty	1 minute	EN 13727 EN 16615*
yeastocidal	with mechanical action dirty	1 minute	EN 13624 EN 16615*
tuberculocidal (M. terrae), mycobactericidal (M. terrae, M. avium)	dirty	1 minute	EN 14348
fungicidal (C. albicans, A. brasiliensis)	dirty	5 minutes	EN 13624 EN 13697
active against enveloped viruses DVV/RKI (Vaccinia strain Elstree, BVDV strain NADL)	10% FBS	1 minute	DVV/RKI
active against enveloped viruses EN 14476 (Modified Vaccinia Ankara)	clean	1 minute	EN 14476
limited spectrum of virucidal activity DVV/RKI (Murine Norovirus, Adenovirus Typ 5)	10% FBS	1 minute	DVV/RKI
virucidal effective against Polyoma SV-40 DVV/RKI	10% FBS	5 minutes	DVV/RKI
virucidal (Murine Norovirus, Adenovirus type 5, Poliovirus type 1)	clean	5 minutes	EN 14476 EN 16777

* tested on reference wipe

Information
about this product



For further information,
please visit our website
www.drdeppe.de

For specific inquiries,
please contact
info@drdeppe.de

CE 0044



Laboratorium Dr. Deppe GmbH
Hooghe Weg 35 · 47906 Kempen
T +49 (0) 2152 5565-0 W www.drdeppe.de

021-06 30.10.2026