Beta Guard rfu wipes

Product information

Product description

Beta Guard rfu wipes are alcohol- and aldehyde-free surface disinfectant wipes in ready-to-use form. The wipe material used is polypropylene.

Area of application

For disinfecting wipe-resistant medical device surfaces, equipment surfaces and objects in the patient environment, as well as surfaces in industrial and food preparation areas. Surfaces of instruments resp. semi-critical and critical medical devices are excluded.



Spectrum of activity

- bactericidal (incl. MRSA), yeasticidal
- tuberculocidal
- active against enveloped viruses (incl. HBV, HIV, HCV, vaccinia-, influenza- and corona-viruses)
- limited spectrum of virucidal activity (incl. noro-, rota-, adeno-viruses)

Instructions for use

The ready-to-use product is used by wiping. It is essential to keep the surface to be disinfected moist throughout the entire exposure times to achieve effectiveness. Surfaces and equipment that come into contact with food should be rinsed thoroughly with drinking-quality water once the exposure time has finished.

Product range and specification

item number	container size	wipe size	container type	packaging unit	PCN (each)	EAN/GTIN (each)	EAN/GTIN (packaging unit)
601111	150 wipes	15 x 18 cm	can	12 x 150 wipes	18312295	4260181253992	4260181251110
601113	150 wipes	15 x 18 cm	can-refill pack	12 x 150 wipes	18312303	4260181254005	4260181251127



DrDeppe

Disinfection wipes

Beta Guard rfu wipes

Product information

Biocide registration number

N-83979
Use disinfectant safely. Always read the label and product information before use.

Composition

The active ingredients contained in 100 g soaking solution:

0.1 g 2-Phenoxyethanol
0.35 g Alkyl(C12-16) dimethyl benzyl ammonium chloride
0.25 g N-(3-Aminopropyl)-N-dodecylpropane-1,3-diamine

Storage and shelf life

Shelf life:
36 months

Shelf life after opening:
28 days

Close the container carefully after each use.
Only use moist wipes.

Store the product in a well-ventilated place in a sealed condition and protect it from unauthorized access.

Exposure times

Beta Guard rfu wipes (polypropylene) spectrum of activity	condition (dirty/ clean)	exposure time	norm
bactericidal	dirty	1 minute	EN 13727* EN 16615
yeasticidal	dirty	1 minute	EN 13624* EN 16615
tuberculocidal (M. terrae)	dirty	60 minutes**	EN 14348*
active against enveloped viruses	dirty	1 minute	EN 14476*
Modified Vaccinia Ankara	clean	1 minute	mod. EN 16615
limited spectrum of virucidal activity (Adenovirus type 5, Murine Norovirus)	dirty	60 minutes**	EN 14476*

* tested on soaking solution
** not for near patient surfaces

Information
about this product

For further information, please visit our website **www.drdeppe.de**

For specific inquiries, please contact **info@drdeppe.de**

DrDeppe

Laboratorium Dr. Deppe GmbH
Hooghe Weg 35 · 47906 Kempen
T +49 (0) 2152 5565-0 **W** www.drdeppe.de

023-05 03.11.2026