Instru Sept AF

Product information

Product description

Instru Sept AF is a concentrated aldehyde- and phenol-free disinfectant for manual disinfection. It is protein-fixing and exhibits very good material compatibility.

Area of application

Instru Sept AF is intended for the reprocessing of medical instruments using the immersion method. Instru Sept AF is not suitable for the final disinfection of semicritical instruments.



Spectrum of activity

- bactericidal (including MRSA), yeasticidal
- tuberculocidal, mycobactericidal
- fungicidal
- active against enveloped viruses (including HBV, HIV, HCV, vaccinia-, influenza-, and corona viruses)
- limited spectrum of virucidal activity (including noro-, rota-, adeno-viruses)

Instructions for use

Instru Sept AF exhibits protein-fixing properties and solubilizes even dried protein residues. It is supplied as a concentrate for preparing a working solution. The recommended concentration for the working solution is 3-5%. Insert instruments after cleaning. All surfaces and cavities must be covered or filled. Please follow the application concentration and contact times. Remove instruments, rinse thoroughly with water, and dry. Follow the hygiene requirements for the reprocessing of medical devices (in Germany: RKI, KRINKO, MPBetreibV). Follow the instructions of the instrument and equipment manufacturers. Do not combine amine-containing Instru Sept AF with aldehyde-based products. Clean and rinse all instruments, accessories, trays, processing equipment, and application aids thoroughly when using the product for the first time or when changing products. For professional use only.

Product range and specification

item number	container size	container type	packaging unit	PCN (each)	EAN/GTIN (each)	EAN/GTIN (packaging unit)
601491	2 L	canister	6 × 2 L	18313337	4260181252537	4260181252544
601492	5 L	canister	2 × 5 L	18313320	4260181252551	4260181252568

Instru Sept AF

Product information

Medical device class

IIb

Composition

The active ingredients contained in 100 g:
8.0 g Didecyldimethylammonium chloride
8.0 g Alkylamines

Storage and shelf life

Shelf life:
36 months

Shelf life after opening:
12 months

Close the container carefully after each use.

Store the product in a well-ventilated place in a sealed condition and protect it from unauthorized access.

Exposure times

Instru Sept AF spectrum of activity	condition (dirty/ clean)	concen-tration	exposure time	norm
bactericidal	dirty	3%	5 min.	EN 13727 EN 14561
yeasticidal	dirty			EN 13624 EN 14562
fungicidal (A. brasiliensis)	dirty	5%	60 min.	EN 13624 EN 14562
tuberculocidal (M. terrae)	dirty	1%	15 min.	EN 14348 EN 14563
mycobactericidal (M. terrae, M. avium)	dirty			EN 14348 EN 14563
active against enveloped viruses (Modified Vaccinia Ankara)	dirty	3%	15 min.	EN 14476 EN 17111
limited spectrum of virucidal activity (Adenovirus type 5, Murine Norovirus)	dirty	5%	30 min.	EN 14476 EN 17111

Information
about this product

For further information, please visit our website **www.drdeppe.de**

For specific inquiries, please contact **info@drdeppe.de**

CE 0044

Laboratorium Dr. Deppe GmbH
Hooghe Weg 35 · 47906 Kempen
T +49 (0) 2152 5565-0 **W** www.drdeppe.de

083-03 09.09.2025