Instru Star rfu

Product information

Product description

Instru Star rfu is a ready-to-use, aldehyde-free instrument disinfectant.

Area of application

It is suitable for the disinfection of medical, surgical instruments using the immersion method. It is also suitable for rotary dental instruments (drills). Not suitable for final disinfection of semicritical instruments.



Spectrum of activity

- bactericidal (including MRSA), yeasticidal
- tuberculocidal, mycobactericidal
- active against enveloped viruses (including HBV, HIV, HCV, vaccinia-, influenza-, and corona viruses)
- virucidal effective against Polyoma SV-40

Instructions for use

Immersion method

Immerse instruments immediately after cleaning in the undiluted liquid. All surfaces and cavities must be covered or filled. Please respect the contact times.

Remove instruments, rinse thoroughly with water, and dry. Do not insert any aluminium parts. Follow the hygiene requirements for the reprocessing of medical devices (in Germany: RKI, KRINKO, MPBetreibV). Follow the instructions of the instrument and equipment manufacturers. For professional use only.

Product range and specification

item number	container size	container type	packaging unit	PCN (each)	EAN/GTIN (each	EAN/GTIN (packaging unit)
601490	2 L	canister	6 × 2 L	18313372	4260181252407	4260181252414
600091	5 L	canister	2 × 5 L	18313366	4260181253947	4260181251363
600071	10 L	canister	1 × 10 L	18313343	4260181253930	4260181251356

p. 2

DrDeppe

Instrument disinfection

Instru Star rfu

Product information

Medical device class

IIb

Composition

The active ingredients contained in 100 g:
1.125 g cocospropylene diammonium borate

Storage and shelf life

Shelf life:
36 months

Shelf life after opening
12 months

Close the container carefully after each use.

Store the product in a well-ventilated place in a sealed condition and protect it from unauthorized access.

Exposure times

Instru Star rfu spectrum of activity	condition (dirty/ clean)	exposure time	norm
bactericidal	dirty	15 min.	EN 13727 EN 14561 (DGHM)
yeasticidal			EN 13624 EN 14562 (DGHM)
tuberculocidal (M. terrae)			EN 14348 EN 14563 (DGHM)
mycobactericidal (M. terrae, M. avium)			EN 14348 EN 14563 (DGHM)
active against enveloped viruses (Vaccinia strain Elstree)			EN 14476
virucidal effective against Polyoma SV40	10% FBS		DVV/RKI

Information
about this product

For further information, please visit our website **www.drdeppe.de**

For specific inquiries, please contact **info@drdeppe.de**

DrDeppe

CE 0044

Laboratorium Dr. Deppe GmbH
Hooghe Weg 35 · 47906 Kempen
T +49 (0) 2152 5565-0 **W** www.drdeppe.de

059-03 15.09.2025