Instru Zym wipes

Product information

Product description

Instru Zym wipes are ready-to-use enzyme-based cleaning wipes. The wipe material is made from cellulose fibres, which are produced from wood, a renewable raw material.

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Instructions for use

After use, wipe the instrument to be cleaned thoroughly with an Instru Zym wipe as soon as possible. If necessary, repeat with a new wipe until the surface is visually clean. Preferably, pre-cleaning with Instru Zym wipes should be carried out from the clean end of the instrument towards the dirty end. For professional use only.

Product range and specification

item number	container size	wipe size	container type	packaging unit	PCN (each)	EAN/GTIN (each)	EAN/GTIN (packaging unit)
601795	80 wipes	20 × 20 cm	flowpack (with LID)	6 × 80 wipes	19691834	4260181257563	4260181257570

Instru Zym wipes

Product information

Instrument pre-cleaning

Medical device class

I

Composition

Ingredients
protease, amylase, lipase, non-ionic surfactants, preservatives, fragrances, purified water

Storage and shelf life

Shelf life:
36 months

Shelf life after opening:
28 days

Close the container carefully after each use.
Only use moist wipes.

Store the product in a well-ventilated place in a sealed condition and protect it from unauthorized access.

Information
about this product

For further information, please visit our website **www.drdeppe.de**

For specific inquiries, please contact **info@drdeppe.de**

CE

Laboratorium Dr. Deppe GmbH
Hooghe Weg 35 · 47906 Kempen
T +49 (0) 2152 5565-0 **W** www.drdeppe.de

073-04 04.11.2026