Labomat MA

Product information

Product description

Labomat MA is a ready-to-use, liquid, mildly alkaline instrument cleaner.

Area of application

Labomat MA is used for machine-based alkaline cleaning.



Application concentration

3 ml / 1 L

Instructions for use

The dosage is automatically applied via the dosing pump of the dishwasher. A temperature of 35-40°C is recommended for the inlet water. Temperatures significantly below 35°C may cause foaming. The instructions of the equipment manufacturers regarding installation, operation, and maintenance must be followed, especially regarding the quality requirements of the process water to be used. After cleaning with Labomat MA, it is recommended to rinse with a neutralizing rinsing agent. For professional use only.

Product range and specification

item number	container size	container type	packaging unit	PCN (each)	EAN/GTIN (each)	EAN/GTIN (packaging unit)
600037	5 L	flat-top canister	3 × 5 L	18312415	4260181251615	4260181251622

DrDeppe

Instrument cleaner

Labomat MA

Product information

Product class

Detergent

Composition

Ingredients:
Chelating agents, polycarboxylates, anionic surfactants, non-ionic surfactants, solubilizers, purified water

Storage and shelf life

Shelf life:
36 months

Shelf life after opening:
12 months

Close the container carefully after each use.

Store the product in a well-ventilated place in a sealed condition and protect it from unauthorized access.

Information
about this product

For further information, please visit our website **www.drdeppe.de**

For specific inquiries, please contact **info@drdeppe.de**

DrDeppe

Laboratorium Dr. Deppe GmbH
Hooghe Weg 35 · 47906 Kempen
T +49 (0) 2152 5565-0 **W** www.drdeppe.de

068-03 08.09.2025