Spray In GA

Product information

Product description

Spray In GA is a ready-to-use, liquid, low-alcoholic disinfectant for surface disinfection. The product shows very good material compatibility under the specified exposure times. Spray In GA is available in the variants Neutral (without fragrance), Melon (fruity fragrance) and Active (fresh fragrance).

Area of application

For rapid disinfection of alcohol resistant surfaces in the vicinity of patients, in clinics, surgeries, out-patient departments and care areas, as well as in industry and public facilities.



Spectrum of activity

- bactericidal (incl. MRSA), yeasticidal
- tuberculocidal, mycobactericidal
- active against enveloped viruses (incl. HBV, HIV, HCV, vaccinia-, influenza- and corona-viruses)
- limited spectrum of virucidal activity (incl. noro-, rota-, adeno-viruses)
- effective against Polyoma SV-40

Instructions for use

The ready-to-use product is preferably to be used in the wiping method using a disposable wipe. Spray disinfection should preferably be limited to areas that are not reachable by wiping disinfection. When using the spray method, it should be combined with wiping method. It is important to keep the surface to be disinfected moist throughout the entire contact time to achieve effectiveness.

Target group (users)

- Nursing staff and care workers
- People in need of care
- Family carers
- Care advisers, independent advice centres or local authority advisers in their advisory capacity

Risks associated with use and contraindications

Spray In GA has no known risks associated with its use or contraindications. There is no evidence of adverse effects or restrictions.

Produktsortiment und Spezifikation

item	container size	container type	packaging unit	PCN (each)	EAN/GTIN (each)	EAN/GTIN (packaging unit)
601939 (Neutral)	500 ml	bottle	18 × 500 ml	20190468	4260181257655	4260181257679
601940* (Melon)	500 ml	bottle	18 × 500 ml	20190474	4260181257693	4260181257716
601941* (Active)	500 ml	bottle	18 × 500 ml	20190445	4260181257730	4260181257754

* on request

p. 2

DrDeppe

Surface disinfection

Spray In GA

Product information

Medical device class

IIa

Composition

The active ingredients contained in 100 g:
20 g Propan-1-ol
0.2 g Alkylbenzyldimethylammonium chloride

Storage and shelf life

Shelf life
36 Monate

Shelf life after opening
12 months

Carefully reclose container after each use.

Store the product tightly closed in a well-ventilate location and protect it against unauthorised access.

Exposure times

Spray In GA spectrum of activity	condition (dirty/ clean)	exposure time	norm
bactericidal	with mechanical action dirty	1 minute	EN 13727 EN 16615*
yeasticidal	with mechanical action dirty	1 minute	EN 13624 EN 16615*
tuberculocidal (M. terrae), mycobactericidal (M. terrae, M. avium)	with mechanical action dirty	2 minutes	VAH 9.1 VAH 14.1*
active against enveloped viruses (Modified Vaccinia Ankara)	with mechanical action clean	1 minute	EN 14476, modified EN16615*
limited spectrum of virucidal activity	clean	30 minutes**	EN 14476
Adenovirus type 5	clean	30 minutes**	EN 14476
Murine Norovirus	clean	1 minute	EN 14476
Rotavirus strain WA	without	30 seconds	EN 14476
virucidal effective against Polyoma SV-40	clean	5 minutes	EN 14476

* tested on reference wipe
** not for near patient surfaces

Informations
about this product

For further information, please visit our website **www.drdeppe.de**

For specific inquiries, please contact **info@drdeppe.de**

CE 0044

DrDeppe

Laboratorium Dr. Deppe GmbH
Hooghe Weg 35 · 47906 Kempen
T +49 (0) 2152 5565-0 **W** www.drdeppe.de

04-0490 02.11.2026