


# Dosing table

for preparing a working solution

Total amount of working solution	Concentration								
	0,25 %	0,5 %	1,0 %	1,5 %	2,0 %	2,5 %	3,0 %	4,0 %	5,0 %
	Quantity of concentrate								
1 litre	2,5 ml	5 ml	10 ml	15 ml	20 ml	25 ml	30 ml	40 ml	50 ml
2 litres	5 ml	10 ml	20 ml	30 ml	40 ml	50 ml	60 ml	80 ml	100 ml
3 litres	7,5 ml	15 ml	30 ml	45 ml	60 ml	75 ml	90 ml	120 ml	150 ml
4 litres	10 ml	20 ml	40 ml	60 ml	80 ml	100 ml	120 ml	160 ml	200 ml
5 litres	12,5 ml	25 ml	50 ml	75 ml	100 ml	125 ml	150 ml	200 ml	250 ml
6 litres	15 ml	30 ml	60 ml	90 ml	120 ml	150 ml	180 ml	240 ml	300 ml
7 litres	17,5 ml	35 ml	70 ml	105 ml	140 ml	175 ml	210 ml	280 ml	350 ml
8 litres	20 ml	40 ml	80 ml	120 ml	160 ml	200 ml	240 ml	320 ml	400 ml
9 litres	22,5 ml	45 ml	90 ml	135 ml	180 ml	225 ml	270 ml	360 ml	450 ml
10 litres	25 ml	50 ml	100 ml	150 ml	200 ml	250 ml	300 ml	400 ml	500 ml

Appilcation example:



Please refer to the label or product information of the relevant product for the required concentration and exposure time. We recommend changing the disinfectant solution daily.



Laboratorium Dr. Deppe GmbH  
 Hooghe Weg 35 · 47906 Kempen  
 T +49 (0) 2152 5565-0 W [www.drdeppe.de](http://www.drdeppe.de)