

Beta Guard rfu Bio-wipes

Product information

Product description

Beta Guard rfu Bio-wipes are ready-to-use, alcohol- and aldehyde-free surface disinfectant wipes. The CO₂ neutral wipe material is 100% biodegradable.

Area of application

For disinfection of wipe-resistant medical device surfaces, equipment surfaces and objects in the vicinity of patients as well as surfaces in the industrial and food sectors. Excluded are surfaces of instruments and semi-critical and critical medical devices.

Spectrum of activity

- bactericidal (incl. MRSA), yeasticidal
- tuberculocidal
- active against enveloped viruses (incl. HBV, HIV, HCV, vaccinia-, influenza- and corona-viruses)
- limited spectrum of virucidal activity (incl. noro-, rota-, adeno-viruses)

Instructions for use

The ready-to-use product is applied by wiping. The surface to be disinfected must be kept moist for the entire exposure time in order to achieve effectiveness.



Product range and specification

item number	container size	wipe size	container type	packaging unit	EAN/GTIN (each)	EAN/GTIN (packaging unit)
601104	80 wipes	20 x 20 cm	flowpack	12 x 80 wipes	4260181254142	4260181251103
601108	70 wipes	29 x 30 cm	bucket	1 x 70 wipes	4260181254128	4260181251080
601109	70 wipes	29 x 30 cm	bucket-refill pack	6 x 70 wipes	4260181254135	4260181251097

Beta Guard rfu Bio-wipes

Product information

Medical device class

IIa

Biocide registration number

N-97387

Use disinfectant safely. Always read the label and product information before use.

Composition

The active ingredients contained in 100 g soaking solution are:

0,1 g 2-Phenoxyethanol

0,35 g Alkyl(C12-16) dimethyl benzyl ammonium chloride

0,25 g N-(3-Aminopropyl)-N-dodecylpropane- 1,3-diamine

Storage and shelf life

Shelf life

flowpack 24 months,

bucket and bucket refill pack 36 months

Shelf life after opening

28 days

Carefully reseal the container after each use.

Only use wet wipes.

Store the product tightly closed in a well-ventilated location and protect it against unauthorised access.

Exposure times

Beta Guard rfu Bio-wipes		condition		norm	exposure time (min.)
		dirty	clean		
bactericidal	with mechanical action	x		EN 13727*	1
				EN 16615	
yeastcidal		x		EN 13624*	
				EN 16615	
tuberculocidal (M. terrae)		x		EN 14348*	60**
active against enveloped viruses		x		EN 14476*	1
Modified Vaccinia Ankara			x	mod. EN 16615	
limited spectrum of virucidal activity (Adenovirus type 5 Murine Norovirus)		x		EN 14476*	60**

* tested on soaking solution

** not for near patient surfaces

Informations

about this product



Further information

will be available on

our homepage

www.drdeppe.de

Special enquiries

should be addressed to

info@drdeppe.de

CE 0044



Laboratorium Dr. Deppe GmbH

Hooghe Weg 35 · 47906 Kempen

T +49 (0) 2152 5565-0 W www.drdeppe.de