

Spray In QF Bio-wipes

Product information

Product description

Spray In QF Bio-wipes are ready-to-use alcoholic surface disinfectant wipes, free from aldehydes and QAV. The CO₂-neutral wipe material is made from renewable raw materials and is 100% biodegradable.

Area of application

For rapid disinfection of alcohol-resistant surfaces in patient environments, clinics, practices, outpatient departments and care areas, as well as in industry and public facilities.

Spectrum of activity

- bactericidal (incl. MRSA), yeasticidal
- tuberculocidal, mycobactericidal
- active against enveloped viruses (incl. HBV, HIV, HCV, vaccinia-, influenza- and corona-viruses)
- limited spectrum of virucidal activity (incl. noro-, rota-, adeno-viruses)

Instructions for use

The ready-to-use product is used by wiping. It is essential to keep the surface to be disinfected moist throughout the entire exposure times to achieve effectiveness.



Product range and specification

item number	container size	wipe size	container type	packaging unit	EAN/GTIN (each)	EAN/GTIN (packaging unit)
601521* **	60 wipes	18 x 20 cm	flowpack	20 x 60 wipes	4260181252926	4260181252933
600971	80 wipes	20 x 20 cm	flowpack	12 x 80 wipes	4260181254050	4260181251202
601102	70 wipes	29 x 30 cm	bucket	1x 70 wipes	4260181254067	4260181251219
601103	70 wipes	29 x 30 cm	bucket-refill pack	6 x 70 wipes	4260181254074	4260181251226

* item is a biocidal product **soaked with 500 ml soaking solution

Spray In QF Bio-wipes

Product information

Medical device class

IIa (except item 601521 – biocidal product only)

Biocide registration number

N-78340

Use disinfectant safely. Always read the label and product information before use.

Composition

The active ingredients contained in 100 g soaking solution:

50 g Ethanol

Storage and shelf life

Shelf life

flowpack 24 months,
bucket and bucket refill pack 36 months

Shelf life after opening

28 days

Close the container carefully after each use.
Only use moist wipes.

Store the product in a well-ventilated place in a sealed condition and protect it from unauthorized access.

Exposure times

Spray In QF Bio-wipes	condition		norm	expo- sure time
	dirty	clean		
bactericidal		x	EN 13727*/ EN 1276* EN 16615	1 min.
yeastical		x	EN 13624*/ EN 1650* EN 16615	
active against enveloped viruses (Modified Vaccinia Ankara)		x	EN 14476* virus- 4-field-test**	
limited spectrum of virucidal activity (Adenovirus type 5, Murine Norovirus)		x	EN 14476*	
tuberculocidal (M. terrae)		x	EN 14348*	5 min.
mycobactericidal (M. terrae und M. avium)	x		EN 14348*	

* tested on soaking solution

** following EN 16615

Informations about this product



For further information,
please visit our website
www.drdeppe.de

For specific inquiries,
please contact
info@drdeppe.de

CE 0044



Laboratorium Dr. Deppe GmbH
Hooghe Weg 35 · 47906 Kempen
T +49 (0) 2152 5565-0 W www.drdeppe.de