

Spray Off Bio-wipes

Product information

Product description

Spray Off Bio-wipes are ready-to-use surface disinfectant wipes that are alcohol- and aldehyde-free. The CO₂-neutral wipe material is made from renewable raw materials and is 100% biodegradable.

Area of application

For surface disinfection of alcohol-resistant surfaces in the patient-near environment, clinics, practices, outpatient clinics and care facilities, as well as in industry and public facilities. Excluded are surfaces of instruments and semi-critical and critical medical devices.

Spectrum of activity

- bactericidal (incl. MRSA), yeasticidal
- tuberculocidal, mycobactericidal
- fungicidal
- active against enveloped viruses (incl. HBV, HIV, HCV, vaccinia-, influenza- and corona-viruses)
- limited spectrum of virucidal activity (incl. noro-, rota-, adeno-viruses)
- virucidal effective against Polyoma SV-40
- virucidal

Instructions for use

The ready-to-use product is used by wiping. It is essential to keep the surface to be disinfected moist throughout the entire exposure times to achieve effectiveness.



Product range and specification

item number	container size	wipe size	container type	packaging unit	EAN/GTIN (each)	EAN/GTIN (packaging unit)
601331	80 wipes	20 x 20 cm	flowpack	12 x 80 wipes	4260181254418	4260181250809
601333	70 wipes	29 x 30 cm	bucket	1 x 70 wipes	4260181254395	4260181250823
601329	70 wipes	29 x 30 cm	bucket refill-pack	6 x 70 wipes	4260181254401	4260181250816

Spray Off Bio-wipes

Product information

Medical device class

IIa

Biocide registration number

N-96077

Use disinfectant safely. Always read the label and product information before use.

Composition

The active ingredients contained in 100 g soaking solution:

63.50 g Ethanol

0.05 g Didecyltrimethylammonium chloride

Storage and shelf life

Shelf life

flowpack 24 months,

bucket and bucket refill-pack 36 months

Shelf life after opening

28 days

Close the container carefully after each use.

Only use moist wipes.

Store the product in a well-ventilated place in a sealed condition and protect it from unauthorized access.

Exposure times

Spray Off Bio-wipes		condition		norm	expo- sure time
		dirty	clean		
bactericidal	with mecha- nical action	x		EN 13727*	1 min.
				EN 16615	
yeastocidal		x		EN 13624*	1 min.
				EN 16615	
tuberculocidal/ mycobactericidal (<i>M. terrae</i> , <i>M. avium</i>)		x		EN 14348*	1 min.
fungicidal (<i>C. albicans</i> , <i>A. brasiliensis</i>)		x		EN 13624*	5 min.
active against enveloped viruses accor- ding to DVV/RKI (Vaccinia strain Elstree, BVDV strain NADL)			10% FBS	DVV/RKI*	1 min.
active against enveloped viruses according to EN 14476 (Modified Vaccinia Ankara)			x	EN 14476	
limited spectrum of viru- cidal activity according to DVV/RKI (Murine Norovirus, Adenovirus type 5)			10% FBS	DVV/RKI*	5 min.
virucidal effective against Polyoma SV-40			10% FBS	DVV/RKI*	
virucidal (Murine Norovirus, Adenovirus type 5, Poliovirus type 1)			x	EN 14476*	

* tested on soaking solution

Informations about this product



For further information,
please visit our website
www.drdeppe.de

For specific inquiries,
please contact
info@drdeppe.de

CE 0044



Laboratorium Dr. Deppe GmbH

Hooghe Weg 35 · 47906 Kempen

T +49 (0) 2152 5565-0 W www.drdeppe.de

