

Lotio 2in1

Product information

Product description

Lotio 2in1 are ready-to-use disinfectant wipes containing Aloe Vera for hand and surface disinfection. The wipe material is made from cellulose fibres, which are produced from the renewable raw material wood.

Area of application

One wipe for two applications: Lotio 2in1 can be used for both, hand disinfection and surface disinfection.



Spectrum of activity

- bactericidal (incl. MRSA), yeasticidal
- active against enveloped viruses (incl. HBV, HIV, HCV, vaccinia-, influenza- and corona-viruses)

Instructions for use

Clean heavily soiled hands and surfaces before disinfection. For hand disinfection, moisten dry hands with the wipe and keep them moist throughout the entire exposure time, if necessary, use another wipe. For surface disinfection, completely wet the surface and if necessary, wipe it again with a second wipe. The surface must be kept moist throughout the entire exposure time. Please mind the exposure times.

Product range and specification

item number	container size	wipe size	container type	packaging unit	EAN/GTIN (each)	EAN/GTIN (packaging unit)
601643	60 wipes	18 x 20 cm	flowpack	20 x 60 wipes	4260181254722	4260181254739


DrDeppe

Lotio 2in1

Product information

Biocide registration number

N-86407

Use disinfectant safely. Always read the label and product information before use.

Composition

The active ingredients contained in 100 g soaking solution:

0.8 g Didecyldimethylammonium chloride

0.5 g 2-Phenoxyethanol

Storage and shelf life

Shelf life

24 months

Shelf life after opening

28 days

Close the container carefully after each use.

Only use moist wipes.

Store the product in a well-ventilated place in a sealed condition and protect it from unauthorized access.

Exposure times

Lotio 2in1	condition		norm	exposure time
	dirty	clean		
bactericidal		x	EN 13727* EN 16615	1 min.
yeastocidal		x	EN 13624* EN 16615	
active against enveloped viruses (Modified Vaccinia Ankara)		x	EN 14476*	

* tested on soaking solution

Informations

about this product



For further information, please visit our website
www.drdeppe.de

For specific inquiries, please contact
info@drdeppe.de



Laboratorium Dr. Deppe GmbH

Hooghe Weg 35 · 47906 Kempen

T +49 (0) 2152 5565-0 W www.drdeppe.de