

Labomat KS

Product information

Product description

Labomat KS is a ready-to-use, phosphate- and surfactant-free rinsing agent for disinfection dishwashers.

Area of application

Labomat KS is used to neutralise alkaline substances on medical instruments and laboratory equipment after fully automated cleaning.

Application concentration

Setting of the dosing pump: 0.1-0.5%

Instructions for use

Dosing takes place automatically via the dosing pump in the dishwasher. Please follow the instructions of the relevant device manufacturer. Do not mix with other products. When changing products, rinse the dosing system including the suction tubing with water. Labomat KS should not be used on objects made of anodised aluminium, chrome-plated or nickel-plated parts. For professional use only.



Product range and specification

item number	container size	container type	packaging unit	PCN (each)	EAN/GTIN (each)	EAN/GTIN (packaging unit)
600021	5 L	flat-top canister	3 × 5 L	18312421	4260181253961	4260181251387

Product class

Detergent

Composition**Ingredients**

Citric acid

Storage and shelf life**Shelf life:**

36 months

Shelf life after opening:

12 months

Close the container carefully after each use.

Store the product in a well-ventilated place in a sealed condition and protect it from unauthorized access.

Information
about this product

For further
information, please
visit our website
www.drdeppe.de

For specific inquiries,
please contact
info@drdeppe.de