

Spray In GA

Product information

Product description

Spray In GA is a ready-to-use, liquid, low-alcoholic disinfectant for surface disinfection. The product shows very good material compatibility under the specified exposure times. Spray In GA is available in the variants Neutral (without fragrance), Melone (fruity fragrance) and Aktiv (fresh fragrance).

Area of application

For rapid disinfection of alcohol resistant surfaces in the vicinity of patients, in clinics, surgeries, outpatient departments and care areas, as well as in industry and public facilities.



Spectrum of activity

- bactericidal (incl. MRSA), yeasticidal
- tuberculocidal, mycobactericidal
- active against enveloped viruses (incl. HBV, HIV, HCV, vaccinia-, influenza- and corona-viruses)
- limited spectrum of virucidal activity (incl. noro-, rota-, adeno-viruses)
- effective against Polyoma SV-40

Instructions for use

The ready-to-use product is preferably to be used in the wiping method using a disposable wipe. Spray disinfection should preferably be limited to areas that are not reachable by wiping disinfection. When using the spray method, it should be combined with wiping method. It is important to keep the surface to be disinfected moist throughout the entire contact time to achieve effectiveness.

Product range and specification

item number	container size	container type	packaging unit	PCN (each)	EAN/GTIN (each)	EAN/GTIN (packaging unit)
601611 (Neutral)	500 ml	bottle	18 × 500 ml	18782859	4260181253459	4260181253466
601612 (Melone)	500 ml	bottle	18 × 500 ml	18782865	4260181253473	4260181253480
601613 (Aktiv)	500 ml	bottle	18 × 500 ml	18782871	4260181253497	4260181253503

Spray In GA

Product information

Biocide registration number

N-109576

Use disinfectant safely. Always read the label and product information before use.

Composition

The active ingredients contained in 100 g:

20 g Propan-1-ol

0.2 g Alkylbenzyltrimethylammonium chloride

Storage and shelf life

Shelf life

36 Monate

Shelf life after opening

12months

Carefully reclose container after each use.

Store the product tightly closed in a well-ventilated location and protect it against unauthorised access.

Informationen about this product



For further information, please visit our website www.drdeppe.de

For specific inquiries, please contact info@drdeppe.de

Exposure times

Spray In GA spectrum of activity	condition (dirty/clean)	exposure time	norm
bactericidal	with mechanical action high	1 minute	EN 13727 EN 16615*
yeasticidal	with mechanical action high	1 minute	EN 13624 EN 16615*
tuberculocidal (M. terrae), mycobactericidal (M. terrae, M. avium)	with mechanical action high	2 minutes	VAH 9.1 VAH 14.1*
active against enveloped viruses (Modified Vaccinia Ankara)	with mechanical action low	1 minute	EN 14476, modified EN16615*
limited spectrum of virucidal activity	low	30 minutes**	EN 14476
Adenovirus type 5	low	30 minutes**	EN 14476
Murine Norovirus	low	1 minute	EN 14476
Rotavirus strain WA	without	30 seconds	EN 14476
virucidal effective against Polyoma SV-40	low	5 minutes	EN 14476

* tested on reference wipe

** not for near patient surfaces

