

Spray In QF wipes

Product information

Product description

Spray In QF wipes are ready-to-use alcohol-based surface disinfectant wipes, free from aldehydes and QAVs. The wipe material is made from cellulose fibres, which are produced from the renewable raw material wood.

Area of application

For rapid disinfection of alcohol-resistant surfaces in patient-related environments, in hospitals, clinics, outpatient clinics and care areas, as well as

in industry and public institutions. This excludes surfaces of instruments resp. semi-critical and critical medical devices. Spray In QF wipes are also suitable for use by persons in the food sector, have been certified as HACCP-compliant in accordance with EU Regulation (EC) No. 852/2004 and can be used safely as part of an HACCP concept.

Spectrum of activity

- bactericidal (incl. MRSA), yeasticidal
- tuberculocidal, mycobactericidal
- active against enveloped viruses (incl. HBV, HIV, HCV, vaccinia-, influenza- and corona-viruses)
- limited spectrum of virucidal activity (incl. noro-, rota-, adeno-viruses)

Instructions for use

The ready-to-use product is used by wiping. It is essential to keep the surface to be disinfected moist throughout the entire exposure times to achieve effectiveness.



Product range and specification

item number	container size	wipe size	container type	packaging unit	PCN (each)	EAN/GTIN (each)	EAN/GTIN (packaging unit)
601799*	60 wipes	18 × 20 cm	flowpack (with LID)	20 × 60 wipes	19692733	4260181257426	4260181257433
601912*	60 wipes	18 × 20 cm	flowpack (with pull tab)	22 × 60 wipes	19881865	4260181250038	4260181250045
601796	80 wipes	20 × 20 cm	flowpack (with LID)	12 × 80 wipes	19692348	4260181257341	4260181257358
601801	70 wipes	29 × 30 cm	bucket	1 × 70 wipes	19692377	4260181257365	4260181257372
601802	70 wipes	29 × 30 cm	bucket-refill pack	6 × 70 wipes	19692578	4260181257389	4260181257396

* item is a biocidal product



Spray In QF wipes

Product information

Medical device class

IIa (except item 601799 and 601912 – biocidal product only)

Biocide registration number

N-115216

Use disinfectant safely. Always read the label and product information before use.

Composition

The active ingredients contained in 100 g soaking solution:

50 g Ethanol

Storage and shelf life

Shelf life:

36 months

Shelf life after opening:

28 days

Close the container carefully after each use.
Only use moist wipes.

Store the product in a well-ventilated place in a sealed condition and protect it from unauthorized access.

Exposure times

Spray In QF wipes (Cellulose) spectrum of activity	condition (dirty/ clean)	exposure time	norm
bactericidal	clean	1 minute	EN 13727* EN 1276* EN 16615
yeasticidal	clean	1 minute	EN 13624* EN 1650* EN 16615
active against enveloped viruses (Modified Vaccinia Ankara)	clean	1 minute	EN 14476*, virus-4-field-test**
limited spectrum of virucidal activity (Adenovirus type 5, Murine Norovirus)	clean	1 minute	EN 14476*
tuberculocidal (M. terrae)	clean	1 minute	EN 14348*
mycobactericidal (M. terrae and M. avium)	dirty	5 minutes	EN 14348*

* tested on soaking solution

** following EN 16615

Information about this product



For further information, please visit our website
www.drdeppe.de

For specific inquiries, please contact
info@drdeppe.de