

Spray Off

Product information

Medical device class

IIa

Biocide registration number

N-91345

Use disinfectant safely. Always read the label and product information before use.

Composition

The active ingredients contained in 100 g:

63.50 g Ethanol

0.05 g Didecyldimethylammonium chloride

Storage and shelf life

Shelf life:

36 months

Shelf life after opening:

12 months

Carefully reclose container after each use.

Store the product tightly closed in a well-ventilated location and protect it against unauthorised access.

Exposure times

Spray Off spectrum of activity	condition (dirty/clean)	exposure time	norm
bactericidal	with mechanical action dirty	1 minute	EN 13727 EN 16615*
yeasticidal	with mechanical action dirty	1 minute	EN 13624 EN 16615*
tuberculocidal (M. terrae), mycobactericidal (M. terrae, M. avium)	dirty	1 minute	EN 14348
fungicidal (C. albicans, A. brasiliensis)	dirty	5 minutes	EN 13624 EN 13697
active against enveloped viruses DVV/RKI (Vaccinia strain Elstree, BVDV strain NADL)	10% FBS	1 minute	DVV/RKI
active against enveloped viruses EN 14476 (Modified Vaccinia Ankara)	clean	1 minute	EN 14476
limited spectrum of virucidal activity DVV/RKI (Murine Norovirus, Adenovirus Typ 5)	10% FBS	1 minute	DVV/RKI
virucidal effective against Polyoma SV-40 DVV/RKI	10% FBS	5 minutes	DVV/RKI
virucidal (Murine Norovirus, Adenovirus type 5, Poliovirus type 1)	clean	5 minutes	EN 14476 EN 16777

* tested on reference wipe

Information
about this product



For further information,
please visit our website
www.drdeppe.de

For specific inquiries,
please contact
info@drdeppe.de

CE 0044



Laboratorium Dr. Deppe GmbH
Hooghe Weg 35 · 47906 Kempen
T +49 (0) 2152 5565-0 W www.drdeppe.de