Instru Suc

Product information

Product description

Instru Suc is a universal, economical, low-foaming and aldehyde-free concentrate for disinfecting suction units.

Area of application

Instru Suc is used for disinfecting suction units with and without amalgam separators.



Spectrum of activity

- bactericidal (including MRSA), yeasticidal

Instructions for use

Biofilm prevention

Daily disinfection with a 2% working solution can effectively prevent biofilm formation.

Biofilm removal

Existing biofilms can also be removed with Instru Suc: After 16 hours of contact time (overnight) with a 2% working solution, no germs were detected for an existing biofilm.

For professional use only.

Product range and specification

item number	container size	container type	packaging unit	PCN (each)	EAN/GTIN (each)	EAN/GTIN (packaging unit)
601489	2 L	canister	6 × 2 L	18313389	4260181252421	4260181252438
600047	5 L	canister	2 × 5 L	18313395	4260181250700	4260181251295

DrDeppe

Instrument disinfection

Instru Suc

Product information

Medical device class

IIb

Composition

The active ingredients contained in 100 g:
7.5 g Dioctyldimethylammoniumchloride

Storage and shelf life

Shelf life:
36 months

Shelf life after opening:
12 months

Close the container carefully after each use.

Store the product in a well-ventilated place in a sealed condition and protect it from unauthorized access.

Exposure times

Instru Suc spectrum of activity	condition (dirty/ clean)	exposure time	concen-tration	norm
bactericidal	dirty	15 min.	2%	EN 13727 (VAH) EN 14561 (VAH 15)
yeasticidal	dirty			EN 13624 (VAH) EN 14562 (VAH 15)

Information
about this product

For further information, please visit our website **www.drdeppe.de**

For specific inquiries, please contact **info@drdeppe.de**

DrDeppe

CE 0044

MADE IN GERMANY 100 %

Laboratorium Dr. Deppe GmbH
Hooghe Weg 35 · 47906 Kempen
T +49 (0) 2152 5565-0 **W** www.drdeppe.de

057-03 15.09.2025